



## PRODUCT SPECIFICATION

catalog number:	<b>D50044</b>
product description:	K2 Polyclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	K2-BgG
source:	Rabbit
form:	Liquid
purification method:	Protein A
buffer:	0.05M phosphate, 0.85% NaCl, pH 7.15
preservative:	0.05% NaN <sub>3</sub>
concentration:	[lot specific] mg/mL
pairing:	This antibody can be paired with K2-BSA V50045314.
notes/applications:	This conjugate is known to work well in lateral flow.
storage:	2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.