



## PRODUCT SPECIFICATION

|                      |   |
|----------------------|---|
| catalog number:      | <b>D50045</b>   |
| product description: | K2 Polyclonal Antibody, Purified  |
| lot number:          | [lot specific]  |
| immunogen:           | K2-BgG  |
| source:              | Rabbit  |
| form:                | Liquid  |
| purification method: | Protein A   |
| buffer:              | 0.05M phosphate, 0.85% NaCl, pH 7.15  |
| preservative:        | 0.05% NaN <sub>3</sub>  |
| concentration:       | [lot specific] mg/mL  |
| pairing:             | This antibody can be paired with K2-BSA V50045314.  |
| notes/applications:  | This conjugate is known to work well in lateral flow.   |
| storage:             | 2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles.   |
| hazard/biohazard:    | This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product. |