



## PRODUCT SPECIFICATION

|                      |   |
|----------------------|---|
| catalog number:      | <b>D51720</b>   |
| product description: | Zolpidem Polyclonal Antibody, Purified  |
| lot number:          | [lot specific]  |
| source:              | Rabbit  |
| form:                | Liquid  |
| purification method: | Protein A purified IgG fraction   |
| buffer:              | 0.05M phosphate, 0.85% NaCl, pH 7.15  |
| preservative:        | none  |
| concentration:       | [lot specific] mg/mL  |
| pairing:             | This antibody can be paired with V51710, zolpidem-HRP or V51720, zolpidem-BSA.  |
| notes/applications:  | This polyclonal is known to work well for lateral flow and ELISA.   |
| storage:             | 2-8°C short term ( $\leq$ 2 weeks); -20°C long term avoiding multiple freeze thaw cycles.   |
| hazard/biohazard:    | This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product. |