



PRODUCT SPECIFICATION

| | |
|----------------------|---|
| catalog number: | D51730 |
| product description: | Zolpidem Polyclonal Antibody, Purified |
| lot number: | [lot specific] |
| source: | Rabbit |
| form: | Liquid |
| purification method: | Protein A purified IgG fraction |
| buffer: | 0.05M phosphate, 0.85% NaCl, pH 7.15 |
| preservative: | 0.05% NaN ₃ |
| concentration: | [lot specific] mg/mL |
| pairing: | This antibody can be paired with V51710, zolpidem-HRP or V51720, zolpidem-BSA. |
| notes/applications: | This polyclonal is known to work well for lateral flow and ELISA. |
| storage: | 2-8°C short term (≤ 2 weeks); -20°C long term avoiding multiple freeze thaw cycles. |
| hazard/biohazard: | This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product. |