



## PRODUCT SPECIFICATION

catalog number:	<b>D53030</b>
product description:	Amphetamine Polyclonal Antibody, Purified
lot number:	[lot specific]
source:	Sheep
form:	Liquid
purification method:	Protein G chromatography
buffer:	0.01M sodium phosphate, 0.15M NaCl (PBS), pH 7.4
preservative:	0.1% NaN <sub>3</sub>
concentration:	[lot specific] mg/mL
expiration date:	[lot specific]
storage:	2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles.
country of origin:	USA
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.