



PRODUCT SPECIFICATION

CATALOG NUMBER: D57010

| | |
|----------------------|--|
| PRODUCT DESCRIPTION: | Rabbit Anti-Phenobarbital Serum |
| HOST ANIMAL: | Rabbit |
| IMMUNOGEN: | Phenobarbital |
| FORM: | Freon treated neat serum |
| PRESERVATIVE: | 0.1% Sodium azide |
| STORAGE TEMPERATURE: | -20°C long term, 2-8°C short term |
| HAZARD/BIOHAZARD: | This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product. |