



## PRODUCT SPECIFICATION

catalog number:	<b>D60230</b>
product description:	<i>Campylobacter jejuni</i> Polyclonal Antibody, Purified
lot number:	[lot specific]
immuogen:	ATCC strain 29428
source:	Rabbit
form:	Purified IgG Fraction
purity:	> 95%
buffer:	0.01M PBS, pH 7.2
preservative:	0.1% NaN <sub>3</sub>
concentration:	[lot specific] mg/mL
titer:	>1:1,000 by indirect ELISA
specificity:	This antibody is specific for whole cells of <i>C. jejuni</i> .
cross-reactivity:	This polyclonal has not been absorbed and may react with related microorganisms.
storage:	2-8°C short term (<6 months); -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.