



PRODUCT SPECIFICATION

CATALOG NUMBER: D60510

PRODUCT DESCRIPTION:	Rabbit Anti <u>Mycobacterium tuberculosis</u> , IgG Fraction
IMMUNOGEN:	Purified PPD
PURITY:	Purified IgG fraction, purity >95%
BUFFER:	0.01M PBS, pH 7.2
PRESERVATIVE:	0.1% NaN ₃
TITER:	neat for gel techniques
STORAGE TEMPERATURE:	2-8°C short term, -20°C long term
SPECIFICITY:	Minimum of 2 major Mycobacterium tuberculosis bands by immunoelectrophoresis (gamma and beta)
OTHER:	This antiserum has not been absorbed and may react with related microorganisms.
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product