



## PRODUCT SPECIFICATION

**CATALOG NUMBER: D60520**

PRODUCT DESCRIPTION:	Poliovirus Polyclonal
HOST ANIMAL:	Goat
IMMUNOGEN:	Pool of purified virus (types 1-3)
FORM:	Neat antisera
PURITY:	>95%
PURIFICATION:	Purified IgG fraction
BUFFER:	0.01M phosphate saline buffer, pH 7.2
PRESERVATIVE:	0.1% NaN <sub>3</sub>
STORAGE TEMPERATURE:	Short term: 2-8°C; Long term: -20°C