



## PRODUCT SPECIFICATION

**CATALOG NUMBER: D60540**

PRODUCT DESCRIPTION:	Anti- <u>Salmonella</u> sp., IgG fraction
HOST ANIMAL:	Rabbit
IMMUNOGEN:	Mixture of <u>S. enteritidis</u> , <u>S. typhimurium</u> , and <u>S. heidelberg</u>
BUFFER:	0.01M PBS, pH 7.2
PRESERVATIVES:	0.1% NaN <sub>3</sub>
TITER:	>1:1,000 by indirect immunofluorescence
STORAGE TEMPERATURE:	-20°C long term; 4-8°C short term
SPECIFICITY:	Polyvalent for <u>Salmonella</u> "O" & "H" antigens