



PRODUCT SPECIFICATION

CATALOG NUMBER: D60550

PRODUCT DESCRIPTION:	Anti-shigella, IgG fraction
HOST ANIMAL:	Rabbit
IMMUNOGEN:	Mixture of <u>S. boydii</u> , <u>S. flexneri</u> , and <u>S. dysenteriae</u>
PURITY:	>95%
BUFFER:	0.01M PBS, pH 7.2
PRESERVATIVES:	0.1% NaN ₃
TITER:	Not determined
STORAGE TEMPERATURE:	-20°C long term; 2-8°C short term
SPECIFICITY:	Polyvalent for Shigella antigens
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.