



PRODUCT SPECIFICATION

CATALOG NUMBER: D60570

PRODUCT DESCRIPTION:	Rabbit anti <u>E. coli</u>
HOST ANIMAL:	Rabbit
IMMUNOGEN:	Mixture of <u>E. coli</u> serotypes
BUFFER:	0.01M phosphate buffered saline, pH 7.2
PRESERVATIVES:	0.1% NaN ₃
TITER:	>1:2,000 by ELISA vs serotypes 0157:H7, 020, 0125, 055, 0111, & K12
STORAGE TEMPERATURE:	-20°C long term, 4°C short term
SPECIFICITY:	Many "O" and "K" antigenic serotypes of <u>Escherichia coli</u>