



PRODUCT SPECIFICATION

CATALOG NUMBER: D60580

PRODUCT DESCRIPTION:	Escherichia coli O157 Polyclonal
HOST ANIMAL:	Goat
IMMUNOGEN:	Mixture of E. coli serotype
FORM:	Neat antisera
PURITY:	>95%
PURIFICATION:	Purified IgG fraction
BUFFER:	0.01M phosphate saline buffer, pH 7.2
PRESERVATIVE:	0.1% NaN ₃
STORAGE TEMPERATURE:	Short term: 2-8°C; Long term: -20°C
CROSS REACTIVITY:	Antiserum is unabsorbed and does crossreact with related Enterobacteriaceae.