



## PRODUCT SPECIFICATION

**CATALOG NUMBER: D60590**

PRODUCT DESCRIPTION:	Staphylococcus Enterotoxin A
HOST ANIMAL:	Rabbit
IMMUNOGEN:	Purified toxin
FORM:	Neat antisera
PURIFICATION:	Purified IgG fraction
BUFFER:	0.01M phosphate saline buffer, pH 7.2
PRESERVATIVE:	0.1% NaN <sub>3</sub>
STORAGE TEMPERATURE:	Short term: 2-8°C; Long term: -20°C
CROSS REACTIVITY:	Non-reactive with other Staph toxins.