



## PRODUCT SPECIFICATION

**CATALOG NUMBER: D60610**

PRODUCT DESCRIPTION:	Rabbit Anti <u>Treponema pallidum</u> , IgG Fraction
IMMUNOGEN:	<u>T. pallidum</u>
PURITY:	Purified IgG fraction, purity >95%
BUFFER:	0.01M PBS, pH 7.2
PRESERVATIVE:	0.1% NaN <sub>3</sub>
TITER:	>1:3,000 by ELISA
SPECIFICITY:	All antigens
STORAGE TEMPERATURE:	2-8°C short term, -20°C long term
OTHER:	This antiserum has not been absorbed and may react with related microorganisms.
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product