



## PRODUCT SPECIFICATION

**CATALOG NUMBER: D60640**

PRODUCT:	Clostridium Difficile Toxin A sp. Polyclonal
HOST ANIMAL:	Rabbit
SOURCE:	Ascites fluid
IMMUNOGEN:	Extract of Cl. perfringens, Cl. sordellii, Cl. novyi, Cl. septicum, Cl. chauvoei
PURIFICATION:	Purified IgG fraction
PURITY:	>95%
BUFFER:	0.01 M PBS, pH 7.2
PRESERVATIVE:	0.1% sodium azide
STORAGE CONDITIONS:	2-8°C short term; -20 C long term
SPECIFICITY:	Reactive with proteins from common Clostridia species.
APPLICATION:	ELISA, fluorescence microscopy, immunoblotting and immunohistochemistry.
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.