



PRODUCT SPECIFICATION

CATALOG NUMBER: D60650

PRODUCT DESCRIPTION:	Goat anti-Cryptosporidium parvum
IMMUNOGEN:	Oocysts purified from bovine feces
PURITY:	Purified IgG fraction, purity >95%
BUFFER:	0.01M PBS, pH 7.2
PRESERVATIVE:	0.1% NaN ₃
TITER:	>1:1,000 by indirect immunofluorescence
STORAGE TEMPERATURE:	-20°C long term; 4-8°C short term
SPECIFICITY:	Intact Oocysts
CROSS REACTIVITY:	Cross-reacts with a large rod-shaped microorganism found in bovine feces
HAZARD/BIOHAZARD:	This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.