



## PRODUCT SPECIFICATION

**CATALOG NUMBER: D60660**

PRODUCT DESCRIPTION:	Giardia Lamblia
HOST ANIMAL:	Goat
IMMUNOGEN:	Purified cysts
FORM:	Neat antisera
PURITY:	>95%
PURIFICATION:	Purified IgG fraction
BUFFER:	0.01M phosphate saline buffer, pH 7.2
PRESERVATIVE:	0.1% NaN <sub>3</sub>
STORAGE TEMPERATURE:	Short term: 2-8°C; Long term: -20°C
CROSS REACTIVITY:	Crossreacts with a rod-shaped microorganism found in rodent feces.