



## PRODUCT SPECIFICATION

**CATALOG NUMBER: D63010**

PRODUCT DESCRIPTION:	Rabbit anti PSA, IgG Fraction
HOST ANIMAL:	Rabbit
IMMUNOGEN:	PSA isolated from pooled human Seminal plasma
FORM:	Liquid
FORM:	DEAE column chromatography purified Immunoglobulin fraction in
BUFFER:	0.015M KPO <sub>4</sub> Buffer, Ph 7.40, 0.85% NaCl
PRESERVATIVE:	0.1% NaN <sub>3</sub>
STORAGE TEMPERATURE:	2-8°C
SUGGESTED USE:	For solid phase coating
SPECIFICITY:	Human Prostate-specific antigen
CROSS REACTIVITY:	No cross reactivity with Human Albumin, AFP, CEA, Ferritin, hCG, PAP, CA 125, CA 15-3, or CA 19-9.