



PRODUCT SPECIFICATION

CATALOG NUMBER: D65040

PRODUCT:	Guinea Pig anti-Porcine Insulin
HOST ANIMAL:	Guinea Pig
IMMUNOGEN:	Porcine Insulin
FORM:	Liquid
PRESERVATIVE:	None
STORAGE TEMPERATURE:	-20°C long term; 2-8°C short term
CROSS REACTIVITY:	Highly reactive to human insulin. Porcine insulin 100%.
HAZARD/BIOHAZARD:	There is no known biohazard associated with this product.