



PRODUCT SPECIFICATION

CATALOG NUMBER: D75000

PRODUCT DESCRIPTION:	Helicobacter Pylori polyclonal antibody, IgG Fraction
HOST ANIMAL:	Rabbit
IMMUNOGEN:	ATCC strain 43504
PURITY:	>95%
BUFFER:	Phosphate Buffered Saline (0.01M, pH 7.2)
PRESERVATIVE:	0.1% sodium azide
TITER:	>1:5,000 by ELISA
STORAGE TEMPERATURE:	2-8°C
SPECIFICITY:	Whole cells lysate
UNINFECTED CELL REACTIVITY:	N/A
CROSS REACTIVITY:	This antiserum has not been absorbed and may react with related microorganisms.
HAZARD/BIOHAZARD:	This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.