



## PRODUCT SPECIFICATION

**CATALOG NUMBER: D80220**

PRODUCT DESCRIPTION:	Goat Anti-human Serum Albumin
HOST:	Goat
IMMUNOGEN:	Purified human albumin protein
BUFFER:	0.01M Phosphate Buffer, 0.075M sodium chloride, pH 7.0, 0.2 micron filtered
PRESERVATIVE:	0.1% Sodium Azide
TOTAL PROTEIN:	7.9 g/dl
TITER:	3.3 mg/ml by Becker method
STORAGE TEMPERATURE:	2-8°C short term, -20°C long-term (avoid repeated freeze thaw cycles)
HAZARD/BIOHAZARD:	This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.