



PRODUCT SPECIFICATION

CATALOG NUMBER: D80320

| | |
|----------------------|---|
| PRODUCT DESCRIPTION: | Goat Antiserum to Human Apolipoprotein A-II |
| SOURCE: | Goat |
| FORM: | Serum, 0.2 micron filtered |
| BUFFER: | 0.05M Tris-HCl, 0.5 M NaCl, pH 7.5 |
| PRESERVATIVE: | 0.09% NaN ₃ |
| TOTAL PROTEIN: | 6.5 g/dL |
| STORAGE TEMPERATURE: | 2-8°C short term, -20°C long term; |
| SPECIFICITY: | Monospecific precipitin reaction, by IEP, with similar electrophoretic mobility as master lot against 2X pooled human serum and neat pooled human plasma. |
| HAZARD/BIOHAZARD: | This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product. |