



PRODUCT SPECIFICATION

CATALOG NUMBER: D80420

PRODUCT DESCRIPTION:	Anti-Human Apolipoprotein A1 Serum
SOURCE:	Goat
IMMUNOGEN:	Purified human Apolipoprotein A1
FORM:	Serum, 0.2 micron filtered
BUFFER:	Tris-Saline, pH 7.4
PRESERVATIVE:	0.1% NaN ₃
BECKER TITER:	1.44 mg/mL
TOTAL PROTEIN:	8.2 g/dL
STORAGE TEMPERATURE:	-20°C long term; 4-8°C short term
SPECIFICITY:	Monospecific by immunoelectrophoresis when tested against pooled human plasma and twice concentrated pooled human serum.
HAZARD/BIOHAZARD:	This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.