



PRODUCT SPECIFICATION

CATALOG NUMBER: D80520

PRODUCT DESCRIPTION:	Anti-Human Apolipoprotein B
HOST ANIMAL:	Goat
IMMUNOGEN:	Purified Human Apolipoprotein B
FORM:	Serum, 0.2 micron filtered
BUFFER:	Tris-saline, pH 7.4
PRESERVATIVE:	0.1% NaN ₃
BECKER TITER:	23.5 mg/ml
TOTAL PROTEIN:	5.4 g/dl
STORAGE TEMPERATURE:	-20°C long term; 4-8°C short term
SPECIFICITY:	Monospecific by IEP with neat antibody versus 2X pooled human serum and neat pooled human plasma.
HAZARD/BIOHAZARD:	This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.