



PRODUCT SPECIFICATION

CATALOG NUMBER: D80620

PRODUCT:	Anti-Human Apolipoprotein Lp(a), IgG Fraction
HOST ANIMAL:	Goat
IMMUNOGEN:	Human Lp(a)
FORM:	Liquid,
BUFFER:	0.2 micron filtered, in 0.05M Tris, 0.25M NaCl, pH 7.51
PRESERVATIVE:	0.09% NaN ₃
TITER:	21.0 mg/ml by Becker Method
TOTAL PROTEIN:	5.38 g/dL by Biuret Method
STORAGE TEMPERATURE:	2-8°C short term, -20°C long term
SPECIFICITY:	Monospecific when tested against normal human serum and pooled human plasma by IEP