



## PRODUCT SPECIFICATION

**CATALOG NUMBER: D80820**

PRODUCT DESCRIPTION:	Anti-Human Apolipoprotein C-III
HOST ANIMAL:	Goat
IMMUNOGEN:	Human Complement C3
BUFFER:	Defibrinated serum, 0.2 um filtered in 0.2 M Tris, 0.5 M NaCl, pH 7.0
PRESERVATIVES:	0.09% sodium azide
TITER:	38.73 OD/mL
STORAGE TEMPERATURE:	4°C short term, -20°C long term
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.