



PRODUCT SPECIFICATION

CATALOG NUMBER: D80920

PRODUCT DESCRIPTION:	Anti-Human Apolipoprotein E Serum
HOST ANIMAL:	Goat
IMMUNOGEN:	Purified Human Recombinant Apolipoprotein E
FORM:	Serum, 0.2 micron filtered
BUFFER:	0.05M Tris HCL, pH 7.5, 0.5 M NaCl
PRESERVATIVE:	0.1% Sodium Azide
TOTAL PROTEIN:	6.8 g/dL
BECKER TITER:	0.61 mg/ml
STORAGE TEMPERATURE:	-20°C long term; 4-8°C short term
SPECIFICITY:	Specific precipitin reaction with similar electrophoretic mobility as master lot against 2X PHS and Apo concentrate.