



PRODUCT SPECIFICATION

CATALOG NUMBER: D81020

PRODUCT DESCRIPTION:	Goat anti-C3
IMMUNOGEN:	Purified Human C3
BUFFER:	0.05M Tris-HCl, 0.5M NaCl, pH 7.5
PRESERVATIVE:	0.1% NaN ₃
TOTAL PROTEIN:	7.9 g/dL
BECKER TITER:	8.3 mg/mL
STORAGE TEMPERATURE:	-20°C long term; 2-8°C short term
PROCESS:	Defibrinated, delipidized and adsorbed by solid phase chromatography as required.
SPECIFICITY:	Monospecific on IEP using undiluted antiserum tested against normal human plasma and twice concentrated normal human serum.