



## PRODUCT SPECIFICATION

**CATALOG NUMBER: D81200**

PRODUCT DESCRIPTION:	Goat anti-Human Ceruloplasmin
BUFFER:	0.05M Tris, 0.5M NaCl, pH 7.5
PRESERVATIVE:	0.09% Sodium Azide
STORAGE TEMPERATURE:	-20°C long term, 2-8°C short term
PROCESS:	Defibrinated, delipidized and adsorbed by solid phase chromatography as required.
SPECIFICITY:	Specific reaction to Ceruloplasmin by IEP with neat antibody versus PHS.
HAZARD/BIOHAZARD:	This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.