



## PRODUCT SPECIFICATION

**CATALOG NUMBER: D82420**

PRODUCT DESCRIPTION:	Anti-Human IgG Fc
HOST ANIMAL:	Goat
FORM:	Serum, IgG fraction
BUFFER:	PBS. pH 7.2
PRESERVATIVE:	0.09% Sodium Azide
ANTIBODY TITER:	12 mg/ml
TOTAL PROTIEN:	50 mg/ml
STORAGE TEMPERATURE:	-20°C long term; 2-8°C short term
HAZARD/BIOHAZARD:	This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.