



## PRODUCT SPECIFICATION

catalog number:	<b>D82510</b>
product description:	Rabbit Anti-Human IgG (H+L), AP
lot number:	[lot specific]
source:	Rabbit
form:	Liquid
purity:	Isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
buffer:	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6
preservative:	none
concentration:	[lot specific] mg/mL
specificity:	Reacts with heavy chains on human IgG and with light chains common to most human immunoglobulins.
cross-reactivity:	This antibody may cross-react with immunoglobulins from other species
expiration date:	[lot specific]
storage:	2-8°C
country of origin:	USA
hazard/biohazard:	There is no known hazard/biohazard associated with this product.