



PRODUCT SPECIFICATION

| | |
|----------------------|---|
| catalog number: | D82520 |
| product description: | Human IgG (H&L) Polyclonal Antibody |
| lot number: | [lot specific] |
| immunogen: | Purified human IgG |
| source: | Goat |
| form: | Liquid; clear to opaque color |
| purification method: | Defibrinated, delipidized and adsorbed by solid phase chromatography; 0.2 μ filtered |
| buffer: | 0.05M Tris, 0.5M NaCl, pH 7.5 |
| preservative: | 0.1% NaN ₃ |
| titer: | [lot specific] mg/mL |
| total protein: | [lot specific] g/dL |
| specificity: | Monospecific by IEP with neat antibody versus 2xPHS and neat PHP |
| manufacture date: | [lot specific] |
| expiration date: | [lot specific] |
| storage: | 2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles. |
| hazard/biohazard: | This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product. |