



PRODUCT SPECIFICATION

CATALOG NUMBER: D82620

PRODUCT DESCRIPTION:	Goat Anti-Human IgM Serum
ANTIGEN:	Purified Human IgM
BUFFER:	0.05M Tris-HCl, 0.5M NaCl, pH 7.5
PRESERVATIVE:	0.1% Sodium Azide
FILTRATION:	0.2 Micron
TOTAL PROTEIN:	7.1 g/dL
BECKER TITER:	19.0 mg/ml
STORAGE TEMPERATURE:	-20°C long term; 2-8°C short term
PROCESS:	Defibrinated, delipidized and adsorbed by solid phase chromatography as required.
SPECIFICITY:	Monospecific by IEP with neat antibody versus 2X pooled human serum and neat pooled human plasma.
HAZARD/BIOHAZARD:	This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.