



## PRODUCT SPECIFICATION

**CATALOG NUMBER: D82720**

PRODUCT DESCRIPTION:	Anti Human Kappa Light Chains (bound & free) serum, IgG fraction
HOST ANIMAL:	Goat
ANTIGEN:	Purified Human Kappa Chains
BUFFER:	PBS (Phosphate Buffered Saline)
PRESERVATIVES:	0.1% NaN <sub>3</sub>
TOTAL PROTEIN:	30 mg/ml
SPECIFIC TITER:	5.0 mg/ml
STORAGE TEMPERATURE:	-20°C long term; 2-8°C short term