



PRODUCT SPECIFICATION

CATALOG NUMBER: D82820

PRODUCT DESCRIPTION:	Anti-Human Lambda Light Chains (free & bound) serum, IgG fraction
HOST ANIMAL:	Goat
ANTIGEN:	Purified Human Lambda Light Chains
BUFFER:	PBS
PRESERVATIVES:	0.1% NaN ₃
TOTAL PROTEIN:	43 mg/ml
SPECIFIC ANTIBODY TITER:	5 mg/ml
STORAGE TEMPERATURE:	-20°C long term; 4-8°C short term