



## PRODUCT SPECIFICATION

**CATALOG NUMBER: D83020**

PRODUCT DESCRIPTION:	Anti-Human Transferrin Serum
HOST ANIMAL:	Goat
ANTIGEN:	Purified Human Transferrin
BUFFER:	Tris-Saline, pH 7.4
PRESERVATIVES:	0.1% NaN <sub>3</sub>
FILTER:	0.2u
TOTAL PROTEIN:	8.0 g/dL (Biuret)
TITER:	8.9 mg/ml (Becker)
STORAGE TEMPERATURE:	2-8°C short term -20°C long term
SPECIFICITY:	Monospecific by IEP when tested against pooled human plasma and twice concentrated pooled human serum.
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.