



## PRODUCT SPECIFICATION

**CATALOG NUMBER: D83220**

PRODUCT DESCRIPTION:	Anti-Human Alpha-1-Antitrypsin serum
HOST ANIMAL:	Goat
ANTIGEN SOURCE:	Purified human Alpha-1-Antitrypsin
BUFFER:	0.05M Tris-HCl pH 7.5 0.5M NaCl
PRESERVATIVES:	0.1% NaN <sub>3</sub>
BECKER TITER:	2.6 mg/ml
TOTAL PROTEIN:	4.5 g/dL
STORAGE TEMPERATURE:	-20°C long term; 2-8°C short term