



## PRODUCT SPECIFICATION

**CATALOG NUMBER: D83320**

PRODUCT DESCRIPTION:	Anti-Human Antithrombin III Polyclonal
HOST ANIMAL:	Goat
ANTIGEN:	Purified Human Antithrombin III
BUFFER:	0.05M Tris-HCl, 0.5M NaCl, pH 7.5
PRESERVATIVES:	0.1% NaN <sub>3</sub>
TOTAL PROTEIN:	3.0 g/dl
BECKER TITER:	1.8 mg/ml
STORAGE TEMPERATURE:	-20°C long term; 4-8°C short term