



## PRODUCT SPECIFICATION

**CATALOG NUMBER: D91100**

PRODUCT DESCRIPTION:	Goat anti Human Prealbumin Polyclonal
SOURCE:	Goat
IMMUNOGEN:	Purified human Prealbumin
BUFFER:	Tris-Saline, pH 7.4
PRESERVATIVE:	0.1% sodium azide
TITER:	2.16 mg/ml (Becker)
TOTAL PROTEIN:	4.8 g/dL (Biuret)
FILTER:	0.2u
STORAGE TEMPERATURE:	2-8°C short term -20°C long term
SPECIFICITY:	Monospecific by immunoelectrophoresis when tested against pooled human plasma and twice concentrated pooled human serum.
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.