



PRODUCT SPECIFICATION

catalog number:	D94090
product description:	<i>H. pylori</i> Polyclonal Antibody, Purified
lot number:	[lot specific]
immunogen:	Total <i>H. pylori</i> sonicate
source:	Rabbit
form:	Liquid
purification method:	Protein A chromatography
buffer:	10mM PBS, pH 7.4
preservative:	0.05% NaN ₃
isotype:	IgG
concentration:	[lot specific] mg/mL
specificity:	This polyclonal antibody is specific for <i>Helicobacter pylori</i> .
expiration date:	[lot specific]
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.