



## PRODUCT SPECIFICATION

catalog number:	<b>G-101-C</b>
product description:	Human IgG (H+L) Polyclonal Antibody, Purified
lot number:	[lot specific]
source:	Goat
form:	Liquid
purification method:	Affinity purified
buffer:	10 mM sodium phosphate, 0.85% (w/v) NaCl, pH 7.2 (PBS)
preservative:	0.05% (w/v) NaN <sub>3</sub>
concentration:	[lot specific] mg/mL
specificity:	This antibody reacts with heavy chains on human IgG and with light chains on all human immunoglobulins. No reactivity is observed to non-immunoglobulin human proteins.
expiration date:	[lot specific]
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.