



PRODUCT SPECIFICATION

CATALOG NUMBER: G-102-C

PRODUCT	Goat anti-human IgG Fc, affinity pure
FORM:	Liquid
PURITY:	Affinity purified
BUFFER:	10 mM sodium phosphate, 0.85% (w/v) sodium chloride, pH 7.2 (PBS)
PRESERVATIVE:	0.05% (w/v) sodium azide
STORAGE:	2-8 degrees C
SPECIFICITY:	This antibody reacts with heavy chains on human IgG, but not with light chains on immunoglobulins. No reactivity is observed to non-immunoglobulin human serum proteins.
NOTE:	Goat antiserum from which this product is derived is obtained from healthy goats of United States origin under the care of a registered veterinarian.