



## PRODUCT SPECIFICATION

**CATALOG NUMBER: G-105-C**

PRODUCT DESCRIPTION:	Goat anti-human IgM, mu chain, Antibody, Affinity Purified
LOT NUMBER:	Reference
FORM:	Liquid
PURITY:	Affinity purified
BUFFER:	10 mM sodium phosphate, 0.85% NaCl (w/v) sodium chloride, pH 7.2 (PBS)
PRESERVATIVE:	0.05% (w/v) sodium azide
CONCENTRATION:	Reference
STORAGE TEMPERATURE:	2-8° C
SPECIFICITY:	This antibody reacts with heavy chains on human IgM, but not with any other immunoglobulins by IEP. No reactivity is observed to non-immunoglobulin human serum proteins.
NOTE:	Goat antiserum from which this product was derived is obtained from healthy goats of United States origin under the care of a registered veterinarian.
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative. Sodium azide may react with lead and copper plumbing to form explosive metal azides. Flush with large amounts water during disposal. There is no known biohazard associated with this product.