



PRODUCT SPECIFICATION

catalog number:	G82320 = G-107-C
product description:	Human IgE Antibody, Affinity Purified
lot number:	[lot specific]
source:	Goat
form:	Liquid
purification method:	Affinity purified
buffer:	10mM sodium phosphate, 0.85% (w/v) NaCl, pH 7.2 (PBS)
preservative:	0.05% NaN ₃
concentration:	[lot specific] mg/mL (based on A ^{1%} _{280nm} = 13.0)
specificity:	This antibody reacts with heavy chains on IgE.
cross-reactivity:	This polyclonal does not react with other immunoglobulins in human serum
expiration date:	[lot specific]
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.