



PRODUCT SPECIFICATION

catalog number:	G-112-C = G81620015
product description:	C-Reactive Protein Polyclonal Antibody, Affinity Purified
lot number:	[lot specific]
source:	Goat
form:	Liquid
purification method:	Affinity purified
buffer:	10mM sodium phosphate, 0.85% (w/v) NaCl, pH 7.2 (PBS)
preservative:	0.05% (w/v) NaN ₃
concentration:	[lot specific] mg/mL
expiration date:	[lot specific]
storage:	2-8°C short term; -20°C long term avoiding multiple freeze thaw cycles.
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.