



PRODUCT SPECIFICATIONS

CATALOG NUMBER: G-118-C

PRODUCT:	Goat anti-human ventricular myosin LC-1, affinity pure
IMMUNOGEN:	Highly purified human ventricular myosin light chain 1
FORM:	Liquid
PURITY:	Affinity purified
BUFFER:	10 mM sodium phosphate, 0.85% (w/v) sodium chloride, pH 7.2 (PBS)
PRESERVATIVE:	0.05% (w/v) sodium azide
STORAGE CONDITIONS:	2-8 degrees C short term -20 degrees C or colder long term
SPECIFICITY:	This antibody reacts with human cardiac myosin light chain 1 by Western blot analysis. It also cross-reacts with human skeletal myosin light chains.
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.