



## PRODUCT SPECIFICATION

**CATALOG NUMBER: G-122-C**

PRODUCT:	Goat anti-myosin N terminal peptide, affinity pure
FORM:	Liquid
PURITY:	Affinity purified
BUFFER:	10 mM sodium phosphate, 0.85% (w/v) sodium chloride, pH 7.2 (PBS)
PRESERVATIVE:	0.05% (w/v) sodium azide
STORAGE:	2-8 degrees C short term -20 degrees C or colder long term
SPECIFICITY:	This antibody reacts only with an N-terminal peptide in the intact LC-1 molecule that substantially reduces cross reactivity to fast skeletal muscle myosin light chains.