



PRODUCT SPECIFICATION

CATALOG NUMBER: G-125-C

| | |
|---------------|--|
| PRODUCT: | Goat anti-human cardiac myoglobin, affinity pure |
| IMMUNOGEN: | Highly purified human cardiac myoglobin |
| FORM: | Liquid |
| PURITY: | Affinity purified. >95% by SDS-PAGE |
| BUFFER: | 10 mM sodium phosphate, 0.85% (w/v) sodium chloride, pH 7.2 (PBS) |
| PRESERVATIVE: | 0.05% (w/v) sodium azide |
| STORAGE: | 2-8 degrees C |
| SPECIFICITY: | This antibody reacts with cardiac and skeletal forms of myoglobin. |