



PRODUCT SPECIFICATION

catalog number:	G-126-C
product description:	Human PSA Polyclonal Antibody, Purified
lot number:	[lot specific]
source:	Goat
form:	Liquid
purification method:	Antigen affinity purified
purity:	> 95% by SDS-PAGE
buffer:	10mM sodium phosphate, 0.85% (w/v) NaCl, pH 7.2 (PBS)
preservative:	0.05% (w/v) NaN ₃
concentration:	[lot specific] mg/mL (based on A ^{1%} _{280nm} = 13.0)
expiration date:	[lot specific]
storage:	2-8°C
hazard/biohazard:	This material contains sodium azide as a preservative; take appropriate precaution. There is no known biohazard associated with this product.